

ML mini abutment system

Screw-retained hybrid dentures

The ML mini abutment system is the appropriate one for correcting the differences in height and angulation in multiple screwed prostheses, achieving passiveness in the superstructure.

Straight mini abutments

SAFETY AND RELIABILITY

A more resistant new screw
(TH.000.M203)
32 N.cm Recommended torque

RESTORATIONS

The system counts with a variety of abutments:

TEMPORARY



COPING



CASTABLE



PROTECTIVE



VARIETY

A wide range of heights and platforms.

CONNECTIONS

Available for all the ML Implant Systems:
SHc, **SHi**, **SHi** CORTICAL and other compatible brands implants.



ML straight mini abutments

Measures

SHe



Mini abutment straight platform 3,45 ht. 2 SHe	PMUE.3520
Mini abutment straight platform 3,45 ht. 4 SHe	PMUE.3540
Mini abutment straight platform 4 ht. 2 SHe	PMUE.4020
Mini abutment straight platform 4 ht. 3 SHe	PMUE.4030
Mini abutment straight platform 4 ht. 4 SHe	PMUE.4040
Mini abutment straight platform 4 ht. 5 SHe	PMUE.4050
Mini abutment straight platform 5 ht. 2 SHe	PMUE.5020
Mini abutment straight platform 5 ht. 4 SHe	PMUE.5040

SHi SHi
CORTICAL




Mini abutment straight platform 4 slim ht. 1 SHi	PMUI.4110
Mini abutment straight platform 4 slim ht. 2 SHi	PMUI.4120
Mini abutment straight platform 4 ht. 1 SHi	PMUI.4010
Mini abutment straight platform 4 ht. 2 SHi	PMUI.4020
Mini abutment straight platform 4 ht. 3 SHi	PMUI.4030
Mini abutment straight platform 5 ht. 1 SHi	PMUI.5010
Mini abutment straight platform 5 ht. 2 SHi	PMUI.5020

SRI




Mini abutment straight platform 3,5 ht. 1 SRI	PMUI.3510.N
Mini abutment straight platform 3,5 ht. 2 SRI	PMUI.3520.N
Mini abutment straight platform 3,5 ht. 3 SRI	PMUI.3530.N
Mini abutment straight platform 4,3 ht. 1 SRI	PMUI.4310.N
Mini abutment straight platform 4,3 ht. 2 SRI	PMUI.4320.N
Mini abutment straight platform 4,3 ht. 3 SRI	PMUI.4330.N
Mini abutment straight platform 5 ht. 1 SRI	PMUI.5010.N
Mini abutment straight platform 5 ht. 3 SRI	PMUI.5030.N



For the placement of the mini abutments, use the hexagon driver for ratchet wrench of 1,75 mm (DCH.170.M). Apply the 32 N.cm torque wrench to adjust the mini abutment to the dental implant (Do not exceed the recommended torque)

IMPORTANT: the hexagon driver for ratchet wrench of 1,75 mm is exclusively designed for the adjustment of the mini abutment components, not for its transference to the mouth cavity.

Angulated mini abutments




VARIETY
Wide range of heights
and platforms

CONNECTIONS

Available for all the ML Implant Systems:
SHe, **SHi**, **SHi** CORTICAL, and other compatible brands implants.

SYSTEM

The ML angulated mini abutments system is simple and easy to be used. It is made up of two parts the mini abutment's base and the cover:



PRESENTATION

The mini abutment comes mounted in a flexible plastic support to facilitate the insertion and placement.

ML angulated 17° mini abutments

Measures

NOTE: The angulated mini abutment comes together with the MP base, MP cover.



SHe

Mini abutment angulated platform 4 ht 2 SHe	PMAE.4020
Mini abutment angulated platform 4 ht. 3 SHe	PMAE.4030
Mini abutment angulated platform 5 ht. 3 SHe	PMAE.5030



SHi

Mini abutment angulated platform 4 ht. 1 SHi	PMAI.4010
Mini abutment angulated platform 4 ht. 2 SHi	PMAI.4020
Mini abutment angulated platform 4 ht. 3 SHi	PMAI.4030
Mini abutment angulated platform 5 ht. 2 SHi	PMAI.5020

ML mini abutments accessories

Measures



Mini abutment Cap 3,45 platform	PPMU.3500
Mini abutment platform Cap 4 platform	PPMU.4000
Mini abutment platform Cap 5 platform	PPMU.5000

NOTE: Screw (TH.000.M203) is included.



Castable mini abutment 3,45 platform	PPMU.3501
Castable mini abutment 4 platform	PPMU.4001
Castable mini abutment 5 platform	PPMU.5001

NOTE: Screw (TH.000.M203) is included.



Temporary mini abutment 3,45 platform	PLMU.3500
Temporary mini abutment 4 platform	PLMU.4000
Temporary mini abutment 5 platform	PLMU.5000

NOTE: Screw (TH.000.M203) is included.



Coping mini abutment 3,45 platform	PTMU.3501
Coping mini abutment 4 platform	PTMU.4001
Coping mini abutment 5 platform	PTMU.5001

NOTE: Screw (TH.050.M207). is included .



Bronze analog mini abutment 3,45 platform	RMU.3502
Titanium analog mini abutment 3,45 platform	RMU.3501
Bronze analog mini abutment 4 platform	RMU.4002
Titanium analog mini abutment 4 platform	RMU.5002
Bronze analog mini abutment 5 platform	RMU.5001
Titanium analog mini abutment 5 platform	RMU.5002